



HL7 Billing Export

APOLLO ADVANCE INTERFACE

HL7 Billing Export Interface Helps Hospitals Reduce Workload & Identify More Billable Items Faster

The HL7 Billing Export Interface automates the billing process, enabling hospitals to reduce staff workload, cut costs, and increase revenue through quicker and more accurate billing of patients. Additional speed and accuracy in billing can be achieved via the interface, allowing hospitals to bill on all clinical events and orders generated by the hospital information system (HIS). **The benefit to hospitals: the convenience and accuracy of a single and comprehensive automated HL7 export file documenting the patient's entire cardiovascular billable activities from start to end of each clinical event.**

Customizable Interface

The HL7 Billing Export Interface transfers data to the hospital billing system. The interface collects data from the Apollo clinical modules and CardioInventory, and creates HL7 P03 export files with the billable data in Financial Transaction 1 (FT1) format. The hospital can customize the LUMEDX HL7 Billing Interface to work with its current billing practices. Configuration allows the site to configure the fields that hold billable items, for example CPT and ICD codes, medications, physicians, and equipment.

Key HL7 Billing Export Interface Benefits:

- Automates daily billing speeding turnaround of patient billing
- Provides single source for all cardiovascular billable data
- Reduces staff data entry through process automation
- Cuts costs, with fewer resources required to run and maintain the billing process
- Eliminates errors associated with duplicate data entry
- Improves workflow process for physicians, clinicians, and administrators
- Increases revenue because patients are billed faster and more accurately
- Generates an immediate ROI

About LUMEDX: With over 350 heart center clients worldwide, LUMEDX is the market leader in fully integrated cardiovascular information systems and the No. 1 independent integrator of cardiology information solutions. LUMEDX offers the most proven, comprehensive package of clinical information tools, cardiovascular products, and services to help medical institutions enhance quality of patient care, reduce costs, streamline workflow, increase patient volume, and grow revenue.

Apollo features the most extensive range of interfaces to clinical systems and HIS—over 60 today—helping to automate the data collection process.

The result: workflow is streamlined saving clients time and resources, and data integrity is maintained.

Flow Diagram for HL7 Billing Export Interface

